

LIFEVEST® WEARABLE DEFIBRILLATOR EMERGENCY PATIENT MANAGEMENT

What alert sounds and voice prompts are being broadcast?

ALERT:

- Device Silent OR Gong Alert (SINGLE TONE)

VOICE:

- None — device silent
- “Contact physician”
- “Treatment has been given, call your doctor”

STATUS:

- Device is monitoring the patient
- Device may be alerting the patient to follow instructions on the screen

Proceed to
First Responder
Instructions
Below

ALERT:

- Siren Alert (TWO TONE)

VOICE:

- “If patient is not responsive, call for help, perform CPR”
- “Device disabled, call ambulance”

STATUS:

- Device cannot detect ECG or the device has delivered the maximum number of treatments

Proceed to
First Responder
Instructions
Below

ALERT:

- Siren Alert (TWO TONE)

VOICE:

- “Press response buttons to delay treatment”
- “Electrical shock possible, **DO NOT TOUCH PATIENT**”
- “Bystanders do not interfere”

STATUS:

- Device has detected a ventricular arrhythmia
- Device is preparing to treat the patient
- Shock likely
- Stop CPR
- Only the patient should press the response buttons (patient consciousness test)
- Do not touch patient
- Allow device to treat the patient

When siren alert stops or “If patient is not responsive, call for help, perform CPR” is broadcast:

Proceed to
First Responder
Instructions
Below

First Responder Instructions

- Proceed with standard evaluation and treatment measures.
- CPR can be performed as long as the device is not broadcasting “**press the response buttons,**” “**electrical shock possible, do not touch patient,**” or “**bystanders do not interfere.**”
- If external defibrillation is available, a decision can be made to remove the LifeVest and monitor/treat the patient with the external equipment.
- To remove the LifeVest, first pull out the battery, then remove the garment from the patient.

24-hour technical support, please call: 800.543.3267